




# Mattress Pressure Mapping System Sensor Element Options

	Light	Standard	High Performance
Sensing points	352 points (11 x 32 matrix)	928 points (16 x 58 matrix)	1,728 points (27 x 64 matrix)
Active area size	55.3 cm x 160 cm 21.8 in x 63 in	46.7 cm x 169.2 cm 18.4 in x 66.6 in	76 cm x 185 cm 30 in x 73 in
Total area size	70.5 cm x 183 cm 27.8 in x 72 in	60.3 cm x 183 cm 23.7 in x 72 in	89.8 cm x 198.8 cm 35.4 in x 78.3
Built-in Wi-fi	Yes	Yes	Yes



**Mattress Sensor**

**Tiedemann Instruments GmbH & Co. KG**

Zur Maximilianshöhe 6  
82467 Garmisch-Partenkirchen  
Germany

Tel.: +49 8821 3068

E-Mail: [info@Tiedemann-Instruments.de](mailto:info@Tiedemann-Instruments.de)

