

FPT-F1 Friction/Peel Tester can be used to test coefficients of static and kinetic friction of plastic films, sheets, foils, paper, cardboard, PP woven bags, fabric( fabric style test), metal-plastic composite strip/belt for communication cable, convey belt and textiles, as well as peel strength test of adhesive laminated products, medical adhesive bandages, release paper and protection films.



### **Professional Technology**

- 2 test modes are available in this instrument including coefficient of friction test and 180° peel test.
- 7 test speeds and 3 test ranges could be selected as required.
- This instrument conforms to multiple standards of ISO, ASTM and GB and the user could select testing methods as required
- Wide range and high-precision of temperature control system to support tests at different temperatures
- This instrument is controlled by computer and micro processing unit which is convenient for user's operation
- The sliding plane and the sled are treated by degaussing and remanence detection which effectively reduce the system errors
- Professional operating software supports superposition analysis of test curves
- Equipped with RS232 port which is convenient to the data transmission and PC connection
- Supports Lystem<sup>TM</sup> Lab Data Sharing System for uniform and systematic data management

#### **Test Standards**

This instrument conforms to various national and international standards:

ISO 8295, ISO 8510-2, ASTM D1894, ASTM D4917, ASTM D3330, TAPPI T816, TAPPI T549, GB 10006, GB/T 2790, GB/T 2791, GB/T 2792

#### **Applications**

	Plastic Films and	Coefficient of static and dynamic friction tests of plastic films and
	Sheets	sheets
	Paper and Paperboard	Coefficient of static and dynamic friction tests of paper and paperboard
	Textiles, Non-woven	Coefficient of static and dynamic friction tests of textiles, non-woven
Basic	c Fabrics and Woven	fabrics and woven bags
<b>Applications</b>	Bags	Tuories und woven ougs
	Rubber	Coefficient of static and dynamic friction tests of rubber products
	Aluminum Foils and	Coefficient of static and dynamic friction tests of aluminum foils, aluminum foil composite films and other metal products
	Aluminum Foil	
	Composite Films	



	Printing Matters	Coefficient of static and dynamic friction tests of printing matters
	High Temperature Test  Wood and Flooring  Photographic Films	Coefficient of static and dynamic friction tests at higher temperature
		than room temperature
		Coefficient of static and dynamic friction tests of wood and flooring
		Coefficient of static and dynamic friction tests of photographic films
Extended	Screen of Mobile	Coefficient of static and dynamic friction tests of mobile screens
<b>Applications</b>	Phone and Leather	against the leather products
	Adhesives	Peeling test of the adhesive products
	Medical Adhesives	Peeling test of medical adhesives e.g. adhesive bandages and plasters
	Protection Films	Peeling test of protective films for mobile phones or computers
	Magnetic Cards	Peeling test of the films and magnetic cards

# **Technical Specifications**

Specifications	FPT-F1	
Capacity Range	0~5 N, 0~10 N, 0~30 N	
Accuracy	0.5% FS	
Mass of Sled	200g (Standard)	
Mass of Sieu	100 g, 500 g, 1000 g, 1814 g, 2000 g are optional	
<b>Test Speed</b>	50 100 150 200 250 300 500 mm/min	
Temperature	Room Temperature ~ 99.9°C	
Power Supply	AC 220V 50Hz	
Instrument	850 mm (L) x 350 mm (W) x 290 mm (H)	
Dimension		
Net Weight	30 kg	

## **Configurations**

Standard	Mainframe, Peeling Grips, Sled of 200g, Sled of 500g, Professional Software and					
Configurations	Communication Cable					
Optional	Test Plate, Sample Cutter, Customized Sled, Customized Sample Grips and Load Cell (with					
	Different Ranges)					

**Please Note:** Labthink is always dedicated to the innovation and improvement of product performance and function. Therefore, technical specifications are subject to change without further notice. Please visit our website at <a href="https://www.labthink.com">www.labthink.com</a> for the latest updates. Labthink reserves the rights of final interpretation and revision.