



Laboratory Equipment

	THK-01	THK-01H
		
Application	Paper, Carton	Paper, Carton and Foils
Test Range (mm)	0-2 (Standard), 0-6 / 0-12 (Optional)	0-2 (Standard)
Resolution (µm)	0,1	0,1
Test Pressure (kPa)	50 ±1 (Paper)	17.5 ±1 (Foils), 50 ±1 (Paper)
Test Area (mm²)	200 (Paper)	50 (Foils), 200 (Paper)
Test Speed	1 - 25/min	1 - 25/min
Automatic Sample Speed	/	0,1 - 99,9 mm/s
Sample Transport	/	0 - 1000 mm
Standards	GB/T 451.3, ASTM D645, ASTM D374, ASTM D1777, TAPPI T411, ISO 4593, ISO 534, ISO 3034, DIN 53105, DIN 53353, JIS K6250, JIS K6328, JIS K6783, JIS Z1702, BS 3983, BS 4817, GB/T 6672, GB/T 6547	ISO 9073-2, ISO 12625-3, ISO 5084, ISO 4593, ISO 534, ISO 3034, GB/T 6672, GB/T 451.3, GB/T 6547, GB/T 6672, GB/T 451.3, GB/T 6547, GB/T 24218.2, ASTM D645, ASTM D374, ASTM D1777, ASTM D6988, ASTM F2251, ASTM D3652, TAPPI T411, TAPPI T411, DIN 53105, DIN 53353, DIN 53370, JIS K6250, JIS K6328, JIS K6783, JIS Z1702, JIS K6250, JIS K6783, JIS Z1702, BS 3983, BS 4817, BS 2782-6, FEFCO No 3, EN 1942